



A Weekly Review (1 of 2)
08/04/2018 - 14/04/2018 (7 Days)
one, patient


This report is the Weekly Review report. This report generates weekly pump and sensor information based on the dates in Date Range A from the Assessment and Progress report. A total of seven(7) days can be displayed on each page.






When using the MiniMed™ 670G system, and an exit out of Auto Mode, which is a feature of the SmartGuard™ technology occurs, the event will be shown as a grey box with a number on top. These grey shaded boxes expand to the duration your MiniMed 670G system operated in Manual Mode. The numbers on top of the grey box refer to the details of the cause for the exit listed in the Exit Reasons Details table, on the bottom right corner of this report.

It is common to have occasional exits on a report if Auto Mode has just been started. If exits occur too often or almost everyday, a review of pump settings may be required.



This portion of the report is the insulin profile. When the pump is operating in Auto Mode, it delivers auto basal. Auto basal can be easily seen by the pink spikes that continuously appear to match your sensor tracing or Continuous Glucose Monitoring [CGM] line. When the pump exits Auto Mode and is operating in Manual Mode, auto basal now becomes a pre-set basal shown as a solid pink line. The pre-set basal is the programmed basal rate set in the pump before training to initiate Auto Mode. This pre-set basal rate is the back up basal rate when Auto Mode exits occur.



The solid black line is the sensor tracing, or the CGM line. Your CGM records up to 288 sensor glucose (SG) values every day and creates this line to show the effect of insulin, carbs, and activity has on the glucose readings.




Review the consistency of the SG readings. Are they always high or above target? Is the glucose frequently below target? The target in CareLink™ is 3.9-10 mmol/L by default and can be modified in CareLink™ | Preferences for report settings. Any duration and frequency below and above target should be addressed during the next clinic visit or follow up phone call. The solid green line inside the target range is a pump target of 6.7mmol/L. This glucose target cannot be adjusted.

Meal periods are also listed below the time to show when meals occur and when SG values change during mealtimes.

Filled in circle - blood glucose (BG) used to calibrate your sensor.

Empty circle - BG entry from the pump or BG meter.


In this report you can also see when an infusion set change occurred. In this example, you can see a small grey arrow above 9am indicating a set change was done Friday.






The circles in purple paired with a purple bar above them are your bolus insulin. Every time there is a bolus for a correction or a meal, a purple circle with a total amount of insulin given and the number of boluses that were given will be shown in parenthesis. For example, during lunch, two(2) boluses totaling 4.4 units of insulin were given. In each bolus, there is an Active Insulin time. This active insulin curve will follow a bolus and extend to the duration of the Active Insulin time set in your pump settings.

The orange boxes below the purple circles are carb entries entered into the pump and used to calculate meal boluses. In this example, two(2) entries were entered during the lunch meal period. These two(2) entries equaled 40 grams of carbs which were used to calculate the total 4.4 units given. Frequent boluses are appropriate as long as all carb entries are entered into the pump so the pump can correctly calculate the appropriate amount of insulin needed for your glucose readings.





Medtronic

Medtronic Australasia Pty Ltd
2 Alma Road
Macquarie Park , NSW 2113
www.medtronic-diabetes.com.au

Safety Information: CareLink™ software

CareLink™ software is intended for use as a tool to help manage diabetes. The purpose of the software is to take information transmitted from insulin pumps, glucose meters and continuous glucose monitoring systems, and turn it into CareLink™ reports. The reports provide information that can be used to identify trends and track daily activities such as carbohydrates consumed, meal times, insulin delivery, and glucose readings. NOTE: CareLink™ report data is intended for use as an adjunct in the management of diabetes only and NOT intended to be relied upon by itself.

For a listing of indications, contraindications, precautions, warnings and potential adverse events please refer to the instructions for Use.

Aooroved 5665012019© 2018 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company.