

# Medtronic



**Medtronic MiniMed**

18000 Devonshire Street  
Northridge, CA 91325 USA  
1 800 646 4633  
+1 818 576 5555  
www.medtronicdiabetes.com

C €0459



# Unomedical

A ConvaTec Company



Unomedical a/s • Aaholmvej 1-3  
Osted • 4320 Lejre • Denmark • www.infusion-set.com

PLACEHOLDER FOR GTIN


PLACEHOLDER FOR PART BARCODE

M995121A022\_1


2021-07

# Medtronic Extended


Infusion Set •




EN	Disconnecting from Site	JA	
DE		SK	
SV		KO	
NL		P2	
DA		HR	
FR		IS	
FI		RO	
IT		ET	
NO		LT	
PT		UK	
ES		BG	
CS		LV	
EL		SR	
HU		MK	
RU		AL	
SL		ID	
PL		MS	
TR		AR	
ZH		HE	




EN	Reconnecting to Site	JA	
DE		SK	
SV		KO	
NL		P2	
DA		HR	
FR		IS	
FI		RO	
IT		ET	
NO		LT	
PT		UK	
ES		BG	
CS		LV	
EL		SR	
HU		MK	
RU		AL	
SL		ID	
PL		MS	
TR		AR	
ZH		HE	



EN	Reservoir Change without Set Change	JA	
DE		SK	
SV		KO	
NL		P2	
DA		HR	
FR		IS	
FI		RO	
IT		ET	
NO		LT	
PT		UK	
ES		BG	
CS		LV	
EL		SR	
HU		MK	
RU		AL	
SL		ID	
PL		MS	
TR		AR	
ZH		HE	



EN	Removing Infusion Set from Body	JA	
DE		SK	
SV		KO	
NL		P2	
DA		HR	
FR		IS	
FI		RO	
IT		ET	
NO		LT	
PT		UK	
ES		BG	
CS		LV	
EL		SR	
HU		MK	
RU		AL	
SL		ID	
PL		MS	
TR		AR	
ZH		HE	



# Medtronic Extended



EN	See page 2	HU		ET	
DE		RU		LT	
SV		SL		UK	
NL		PL		BG	
DA		TR		LV	
FR		ZH		SR	
FI		JA		MK	
IT		SK		AL	
NO		KO		ID	
PT		P2		MS	
ES		HR		AR	
CS		IS		HE	
EL		RO			

## Indications for use

The Medtronic Extended infusion set is indicated for subcutaneous infusion of medication administered by an external pump. The infusion set is indicated for single use.

## Description

The Medtronic Extended infusion set has a 90-degree soft cannula. It is delivered ready to use in a pre-loaded insertion device with automatic needle retraction. It is provided sterile and nonpyrogenic.

## Intended use

The infusion set is indicated for subcutaneous infusion of insulin in the treatment of diabetes mellitus. Refer to the Medtronic insulin pump user guide for a list of compatible insulins and infusion sets.

## Contraindications

This infusion set is indicated for subcutaneous use only. Do not use the infusion set for intravenous infusion. Do not use the infusion set with blood or blood products.

## General instructions

- The infusion set may be worn for a maximum of seven days, or according to the insulin labeling, whichever duration is shorter.
- If the infusion set is being used for the first time, do the first set-up in the presence of a healthcare professional.

## Warnings

- Do not use the infusion set if the package is opened or damaged. An opened or damaged package may be contaminated. This can cause

an infection.

- Do not use the infusion set beyond seven days as this may cause inaccurate insulin delivery, infection, or site irritation.
- Do not use the infusion set beyond the Use-by date, or if the package is opened or damaged, as sterility may be compromised.
- Do not use the infusion set if the disconnect cover is removed. Use a new infusion set instead. One purpose of the disconnect cover is as a safety feature against accidental firing of the insertion device which could cause injury.
- Do not change the infusion set just before bedtime unless blood glucose can be checked one to three hours afterwards.
- Do not put alcohol, disinfectants, perfumes, deodorants, cosmetics, or other substances with solvents on the infusion set. These substances may damage the integrity of the infusion set.
- Never point the insertion device toward any body part where insertion is not desired.
- Always rotate the insertion site when changing the infusion set. Reusing an insertion site too often may cause scarring and unpredictable insulin delivery. Review the pump user guide on how to rotate insertion sites.
- Check the insertion site often through the clear window. Improper insertion and maintenance of the insertion site can cause inaccurate insulin delivery, infection, or site irritation. Replace the infusion set at a new site if the soft cannula is not properly inserted.
- Check the cannula housing and tubing for blood. Blood may cause insufficient insulin delivery. This may result in high blood glucose. Replace the infusion set at a new site if there is blood.

- Fill the tubing completely with insulin before insertion. Do not leave air in the tubing. An uncontrolled amount of insulin may be delivered if air is present.
- Check the tubing for clogs or leaks if a blood glucose reading is high. Clogs or leaks may restrict insulin delivery and result in high blood glucose. Replace the infusion set even if a clog or leak is suspected, but not found.
- Do not attempt to clear air or an occlusion in the tubing while it is connected to the body. An uncontrolled amount of insulin may be delivered. This can cause high or low blood glucose. Disconnect the tubing before adjusting it.
- Do not reuse the infusion set. Reuse of the infusion set may cause damage to the cannula or needle and may lead to infection, site irritation, or unpredictable insulin delivery.
- If insulin or any other liquid gets inside the tubing connector, it can temporarily block the vents that allow the pump to properly fill the infusion set. **This may result in the delivery of too little or too much insulin, which can cause hyperglycemia or hypoglycemia.** If this occurs, start over with a new reservoir and infusion set.
- For no delivery or Insulin Flow Blocked alarms, refer to the pump user guide.
- Do not use the reservoir or infusion set if either is damaged while connecting. A damaged reservoir or infusion set may cause too little delivery of insulin due to leakage, which may cause hyperglycemia.

## Precautions

- Do not use another type of infusion set without consulting a healthcare professional for correct

handling. A healthcare professional should always be consulted when choosing an infusion set.

- Choose insertion sites as recommended by a healthcare professional. The choice of site depends on treatment and patient-specific factors such as body composition and physical activity level.
- Check that the insertion site is free of skin irritation, such as redness, scar tissue, or bleeding. Do not insert the infusion set into muscle or over bone. Doing so can cause pain or may damage the infusion set. If the infusion set is improperly inserted, replace the infusion set and insert at a new site.
- If needed, remove body hair around the insertion site to ensure that the adhesive sticks to the skin.
- Do not use an insertion site located under a belt or waistband, or where the site is constrained by clothing or accessories. Insertion at these locations can cause disconnection of the infusion set and interruption of insulin delivery, and lead to hyperglycemia.
- Do not reposition the infusion set on the body after the adhesive is placed on the skin. Repositioning the infusion set may damage the adhesive. Replace the infusion set if the adhesive is damaged.
- Always check blood glucose one to three hours after inserting a new infusion set, or after changing the reservoir while using the same infusion set. This will confirm that the amount of insulin being infused is having the intended effect on blood glucose levels.

- Check blood glucose several times throughout the day or as recommended by a healthcare professional.
- Consult a healthcare professional for how to correct for missed insulin and for how much time the pump can remain disconnected.

### Storage and disposal

- Store infusion sets in a cool, dry location at room temperature. Do not store infusion sets in direct sunlight, or in high humidity.
- Store and handle insulin according to the manufacturer's instructions.
- Dispose of the insertion device in a proper sharps container according to local laws.
- Dispose of a used infusion set according to local regulations for biohazardous waste.

### Serious incident

- If, during the use of this device or as a result of its use a serious incident has occurred, please report it to the manufacturer and a national authority.

### Warranty

- For product warranty information, please contact a local Medtronic support representative, or visit:  
[www.medtronicdiabetes.com/warranty](http://www.medtronicdiabetes.com/warranty)

### Additional information for the healthcare professional


- Consider treatment and patient-specific factors when choosing an insertion site.
- The length of the soft cannula should be based on treatment and patient factors such as physiology and activity level. If the soft cannula is too long, insertion pain or insertion into the

bone or muscle may occur. If the soft cannula is too short, leakage or irritation at the insertion site may occur. In both cases, the infusion must be stopped and a new soft cannula with a different length must be inserted at a new insertion site.

- Inform the patient that the selection of an insertion site is dependent on cannula length.

© 2020 Medtronic. All rights reserved. Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company.

## Components



- |  |                                 |
|--|---------------------------------|
| A. Top button                                | H. Paper backing                |
| B. Disconnect cover                          | I. Adhesive                     |
| C. Insertion device                          | J. Tubing connector             |
| D. Raised arrows indicating tubing direction | K. Tubing                       |
| E. Clear window                              | L. Site connector               |
| F. Cannula housing                           | M. White cap for site connector |
| G. Soft cannula                              |                                 |

## Instructions for use

- Read all instructions carefully before using the Medtronic Extended infusion set.
- Consult the pump user guide for important information about pump therapy prior to connecting the infusion set. This information includes connection and filling procedures, possible errors, and potential risks related to pump therapy.
- Follow good hygienic procedures. If this infusion set is being used for the first time, perform the first set-up in the presence of a healthcare

professional.

## Inserting infusion set into body

Unfold the front cover of this booklet to view images that are keyed to the steps below.

1. Wash hands with soap and water.
2. Select a recommended insertion site (shown in gray) and as indicated by a healthcare professional.
3. Clean the insertion site with a disinfectant as directed by a healthcare professional. Allow to air dry before inserting the infusion set. Remove hair around the insertion site to ensure the adhesive sticks to the skin if needed.
4. Open the package and remove the paper from the tubing.
5. Place the tubing connector on the top of a newly filled reservoir. If a MiniMed tubing connector is in use, twist the tubing connector clockwise until it locks in place.

For the MiniMed tubing connector: If insulin or any other liquid gets inside the tubing connector, it can temporarily block the vents that allow the pump to properly fill the infusion set. **This may result in the delivery of too little or too much insulin, which can cause hyperglycemia or hypoglycemia.** If this occurs, start over with a new reservoir and infusion set.

6. Remove the white cap from the site connector by (a) squeezing the sides of the site connector then (b) pulling the white cap off. Keep the white cap for later use.
7. Put the reservoir into the pump, then

completely fill the tubing with insulin as instructed by the pump user guide. Do not leave any air in the tubing. The tubing is completely filled when drops of insulin are seen coming out of the site connector.

8. Remove the paper backing from the adhesive. Be careful not to touch the adhesive.
9. Remove the disconnect cover from the insertion device by gently squeezing the sides of the disconnect cover then pulling it away from the insertion device. Keep the disconnect cover for later use. The raised arrows indicate the tubing direction.
- 10a. Stretch the skin until smooth. Then press the insertion device against the skin.
- 10b. Press the top button down completely to insert the infusion set.
11. Gently and carefully remove the insertion device.
12. Gently press the adhesive onto the skin with one finger. Replace the infusion set if the adhesive does not stick to the skin.
13. Gently hold the cannula housing steady with one finger. Then push the site connector straight into the cannula housing until it clicks.
14. Fill the soft cannula with insulin:

6 mm: 0.6 units (0.006 ml)

9 mm: 0.6 units (0.006 ml)

Dispose of the insertion device in an appropriate sharps container according to local laws. Keep the disconnect cover and the white cap for use when the infusion set is disconnected.

## Disconnecting from site

1. Wash hands with soap and water.
2. Gently hold the cannula housing steady with one finger. Then (a) squeeze the sides of the site connector and (b) pull the site connector out of the cannula housing.
3. Put the white cap on the site connector.
4. Put the disconnect cover on the cannula housing. Then push the cover into the cannula housing until it clicks.

## Reconnecting to site

1. Wash hands with soap and water.
2. Use one finger to gently hold the cannula housing steady. Then (a) squeeze the sides of the disconnect cover and (b) pull the disconnect cover out from the cannula housing.
3. Remove the white cap from the site connector by (a) squeezing the sides of the site connector then (b) pulling the white cap off. Ensure that there is no air in the tubing. **If there are air bubbles in the tubing:** Fill the tubing with insulin as instructed in the pump user guide. Do not leave any air in the tubing.
4. Use one finger to gently hold the cannula housing steady. Then, push the site connector into the cannula housing until it clicks.

## Reservoir change without set change

1. Wash hands with soap and water.
2. Use one finger to gently hold the cannula housing steady. Then (a) squeeze the sides of

the site connector and (b) pull it out from the cannula housing.

3. Remove the reservoir from the pump then (a) remove the tubing connector from the reservoir by turning the connector counter-clockwise. Dispose of the used reservoir in an appropriate container according to local laws. Do not throw away the tubing. Set it aside on a clean surface. For specific instructions on removing reservoir from the pump, refer to the pump user guide. (b) Place the tubing connector on the top of a newly filled reservoir. Twist the tubing connector clockwise until it locks in place.

**For the tubing connector, make sure there is no liquid inside the tubing connector or on top of the reservoir. Liquid can block vents and cause inaccurate insulin flow. This can result in high or low blood glucose. If there is liquid, use a new infusion set and a new reservoir.**

4. Put the reservoir into the pump. Ensure that there is no air in the tubing. **If there are air bubbles in the tubing:** Fill the tubing with insulin as instructed in the pump user guide. Do not leave any air in the tubing.
5. Use one finger to gently hold the cannula housing steady. Then, push the site connector into the cannula housing until you hear a click. **Do not fill the cannula. When reconnecting the infusion set following a reservoir-only change, filling the cannula may result in the delivery of too much insulin, which could cause hypoglycemia.**

## Removing infusion set from body

1. Wash hands with soap and water.
2. Use one finger to gently hold the cannula housing steady. Then (a) squeeze the sides of the site connector and (b) pull the site connector out of the cannula housing.
3. Remove the reservoir from the pump. Remove the tubing connector from the reservoir by turning the connector counter-clockwise. For specific instructions on removing the reservoir from the pump, refer to the pump user guide.
4. Carefully lift the adhesive around the cannula housing. Then, pull the soft cannula out of the skin.









EN		Caution	JA
DE			SK
SV			KO
NL			P2
DA			HR
FR			IS
FI			RO
IT			ET
NO			LT
PT			UK
ES			BG
CS			LV
EL			SR
HU			MK
RU			AL
SL			ID
PL			MS
TR			AR
ZH			HE





EN		Do not re-use	JA
DE			SK
SV			KO
NL			P2
DA			HR
FR			IS
FI			RO
IT			ET
NO			LT
PT			UK
ES			BG
CS			LV
EL			SR
HU			MK
RU			AL
SL			ID
PL			MS
TR			AR
ZH			HE



EN		Manufacturer	JA
DE			SK
SV			KO
NL			P2
DA			HR
FR			IS
FI			RO
IT			ET
NO			LT
PT			UK
ES			BG
CS			LV
EL			SR
HU			MK
RU			AL
SL			ID
PL			MS
TR			AR
ZH			HE

<div style="border: 2px solid black; padding: 2px; display: inline-block;"> <b>LOT</b> </div>			
EN	Batch code	JA	
DE		SK	
SV		KO	
NL		P2	
DA		HR	
FR		IS	
FI		RO	
IT		ET	
NO		LT	
PT		UK	
ES		BG	
CS		LV	
EL		SR	
HU		MK	
RU		AL	
SL		ID	
PL		MS	
TR		AR	
ZH		HE	

			
EN	Keep dry	JA	
DE		SK	
SV		KO	
NL		P2	
DA		HR	
FR		IS	
FI		RO	
IT		ET	
NO		LT	
PT		UK	
ES		BG	
CS		LV	
EL		SR	
HU		MK	
RU		AL	
SL		ID	
PL		MS	
TR		AR	
ZH		HE	

			
EN	Open here	JA	
DE		SK	
SV		KO	
NL		P2	
DA		HR	
FR		IS	
FI		RO	
IT		ET	
NO		LT	
PT		UK	
ES		BG	
CS		LV	
EL		SR	
HU		MK	
RU		AL	
SL		ID	
PL		MS	
TR		AR	
ZH		HE	



EN	Distributed by	JA	
DE		SK	
SV		KO	
NL		P2	
DA		HR	
FR		IS	
FI		RO	
IT		ET	
NO		LT	
PT		UK	
ES		BG	
CS		LV	
EL		SR	
HU		MK	
RU		AL	
SL		ID	
PL		MS	
TR		AR	
ZH		HE	

**MD**

EN	Medical device	JA	
DE		SK	
SV		KO	
NL		P2	
DA		HR	
FR		IS	
FI		RO	
IT		ET	
NO		LT	
PT		UK	
ES		BG	
CS		LV	
EL		SR	
HU		MK	
RU		AL	
SL		ID	
PL		MS	
TR		AR	
ZH		HE	



EN	Non-pyrogenic	JA	
DE		SK	
SV		KO	
NL		P2	
DA		HR	
FR		IS	
FI		RO	
IT		ET	
NO		LT	
PT		UK	
ES		BG	
CS		LV	
EL		SR	
HU		MK	
RU		AL	
SL		ID	
PL		MS	
TR		AR	
ZH		HE	



EN	Consult instructions for use	JA	
DE		SK	
SV		KO	
NL		P2	
DA		HR	
FR		IS	
FI		RO	
IT		ET	
NO		LT	
PT		UK	
ES		BG	
CS		LV	
EL		SR	
HU		MK	
RU		AL	
SL		ID	
PL		MS	
TR		AR	
ZH		HE	



EN	Single sterile barrier system	JA	
DE		SK	
SV		KO	
NL		P2	
DA		HR	
FR		IS	
FI		RO	
IT		ET	
NO		LT	
PT		UK	
ES		BG	
CS		LV	
EL		SR	
HU		MK	
RU		AL	
SL		ID	
PL		MS	
TR		AR	
ZH		HE	



EN	Replace every X days	JA	
DE		SK	
SV		KO	
NL		P2	
DA		HR	
FR		IS	
FI		RO	
IT		ET	
NO		LT	
PT		UK	
ES		BG	
CS		LV	
EL		SR	
HU		MK	
RU		AL	
SL		ID	
PL		MS	
TR		AR	
ZH		HE	


## REF

EN	Catalogue number	JA
DE		SK
SV		KO
NL		P2
DA		HR
FR		IS
FI		RO
IT		ET
NO		LT
PT		UK
ES		BG
CS		LV
EL		SR
HU		MK
RU		AL
SL		ID
PL		MS
TR		AR
ZH		HE



EN	Use-by-date (YYYY-MM-DD)	JA
DE		SK
SV		KO
NL		P2
DA		HR
FR		IS
FI		RO
IT		ET
NO		LT
PT		UK
ES		BG
CS		LV
EL		SR
HU		MK
RU		AL
SL		ID
PL		MS
TR		AR
ZH		HE

GTIN		
EN	Global Trade Item Number	JA
DE		SK
SV		KO
NL		P2
DA		HR
FR		IS
FI		RO
IT		ET
NO		LT
PT		UK
ES		BG
CS		LV
EL		SR
HU		MK
RU		AL
SL		ID
PL		MS
TR		AR
ZH		HE

		
EN	Keep away from sunlight	JA
DE		SK
SV		KO
NL		P2
DA		HR
FR		IS
FI		RO
IT		ET
NO		LT
PT		UK
ES		BG
CS		LV
EL		SR
HU		MK
RU		AL
SL		ID
PL		MS
TR		AR
ZH		HE



EN	Do not use if package is damaged	JA	
DE		SK	
SV		KO	
NL		P2	
DA		HR	
FR		IS	
FI		RO	
IT		ET	
NO		LT	
PT		UK	
ES		BG	
CS		LV	
EL		SR	
HU		MK	
RU		AL	
SL		ID	
PL		MS	
TR		AR	
ZH		HE	

**R<sub>x</sub> Only**

EN	Requires a prescription in the USA	JA	
DE		SK	
SV		KO	
NL		P2	
DA		HR	
FR		IS	
FI		RO	
IT		ET	
NO		LT	
PT		UK	
ES		BG	
CS		LV	
EL		SR	
HU		MK	
RU		AL	
SL		ID	
PL		MS	
TR		AR	
ZH		HE	



**STERILE EO**

EN	Sterilized using ethylene oxide	JA	
DE		SK	
SV		KO	
NL		P2	
DA		HR	
FR		IS	
FI		RO	
IT		ET	
NO		LT	
PT		UK	
ES		BG	
CS		LV	
EL		SR	
HU		MK	
RU		AL	
SL		ID	
PL		MS	
TR		AR	
ZH		HE	



EN	Date of manufacture (YYYY-MM-DD)	JA	
DE		SK	
SV		KO	
NL		P2	
DA		HR	
FR		IS	
FI		RO	
IT		ET	
NO		LT	
PT		UK	
ES		BG	
CS		LV	
EL		SR	
HU		MK	
RU		AL	
SL		ID	
PL		MS	
TR		AR	
ZH		HE	

# C €0459

EN	Conformité Européenne (European Conformity) This symbol means that the device fully complies with applicable European Union Acts.
DE	
SV	
NL	
DA	
FR	
FI	
IT	
NO	
PT	
ES	
CS	
EL	
HU	
RU	
SL	
PL	
TR	
ZH	

# C €0459

JA	
SK	
KO	
P2	
HR	
IS	
RO	
ET	
LT	
UK	
BG	
LV	
SR	
MK	
AL	
ID	
MS	
AR	
HE	

**Contacts:**

**Africa:**  
Medtronic South Africa and Southern Africa  
Office Reception Tel:  
+27(0) 11 260 9300  
Diabetes: 24/7 Helpline:  
0800 633 7867  
Sub-Sahara 24/7 Helpline:  
+27(0) 11 260 9490

**Africa:**  
Medtronic South Africa and Southern Africa  
Office Reception Tel:  
+27(0) 11 260 9300  
Diabetes: 24/7 Helpline:  
0800 633 7867  
Sub-Sahara 24/7 Helpline:  
+27(0) 11 260 9490

**Albania:**  
Net Electronics Albania  
Tel.: +355 697070121

**Argentina:**  
Corpomedica S.A.  
Tel.: +(11) 4 814 1333  
Medtronic Directo 24/7:  
+0800 333 0752

**Armenia:**  
Exiol LLC  
Tel.: +374 98 92 00 11  
or +374 94 38 38 52

**Australia:**  
Medtronic Australasia Pty. Ltd.  
Tel.: 1800 668 670

**Bangladesh:**  
Sonargaon Healthcare Pvt Ltd.  
Mobile: +(91)-9903995417  
or +(880)-1714217131

**Belarus:**  
Zarga Medica  
Tel.: +37517 336 97 00  
+37529 613 08 08  
+37517 215 02 89  
Helpline: +74995830400

**België/Belgique:**  
N.V. Medtronic Belgium S.A.  
Tel.: 0800-90805

**Bosnia and Herzegovina:**  
"Novopharm" d.o.o. Sarajevo  
Tel.: +387 33 476 444  
Helpline: 0800 222 33  
Epsilon Research Intern. d.o.o.  
Tel.: +387 51 251 037  
Helpline: 0800 222 33

**Brasil:**  
Medtronic Comercial Ltda.  
Tel.: +(11) 2182-9200  
Medtronic Directo 24/7:  
+0800 773 9200

**Bulgaria:**  
RSR EOOD  
Tel.: +359 888993083  
Helpline: +359 884504344

**Canada:**  
Medtronic Canada ULC  
Tel.: 1-800-284-4416 (toll free/sans frais)

**Chile:**  
Medtronic Chile  
Tel.: +(9) 66 29 7126  
Medtronic Directo 24/7:  
+1 230 020 9750  
Medtronic Directo 24/7 (From Santiago):  
+(2) 595 2942

**China:**  
Medtronic (Shanghai) Management Co., Ltd.  
Landline: +86 800-820-1981  
Mobile Phone: +86 400-820-1981  
Calling from outside China: +86 400-820-1981

**Colombia:**  
Medtronic Latin America Inc. Sucursal Colombia  
Tel.: +(1) 742 7300  
Medtronic Directo 24/7 (Landline):  
+01 800 710 2170  
Medtronic Directo 24/7  
(Cellular): +1 381 4902

**Croatia:**  
Mediligo d.o.o.  
Tel.: +385 1 6454 295  
Helpline: +385 1 4881144  
Medtronic Adriatic d.o.o.  
Helpline: +385 1 4881120

**Česká republika:**  
Medtronic Czechia s.r.o.  
Tel.: +420 233 059 111  
Non-Stop Helpline (24/7):  
+420 233 059 059  
Zákaznický servis (8:00 - 17:00):  
+420 233 059 950

**Danmark:**  
Medtronic Danmark A/S  
Tel.: +45 32 48 18 00

**Deutschland:**  
Medtronic GmbH  
Geschäftsbereich Diabetes  
Telefon: +49 2159 8149-370  
Telefax: +49 2159 8149-110  
24-Stdn-Hotline: 0800 6464633

**Eire:**  
Accu-Science Ltd.  
Tel.: +353 45 433000

**España:**  
Medtronic Ibérica S.A.  
Tel.: +34 91 625 05 42  
Fax: +34 91 625 03 90  
24 horas: +34 900 120 330

**Estonia:**  
AB Medical Group Estonia Ltd  
Tel.: +372 6552310  
Helpline: +372 5140694

**Europe:**  
Medtronic Europe S.A. Europe,  
Middle East and Africa HQ  
Tel.: +41 (0) 21-802-7000

**France:**  
Medtronic France S.A.S.  
Tel.: +33 (0) 1 55 38 17 00

**Hellas:**  
Medtronic Hellas S.A.  
Tel.: +30 210677-9099

**Hong Kong:**  
Medtronic Hong Kong Medical Ltd.  
Tel.: +852 2919-1300  
To order supplies: +852 2919-1322  
24-hour helpline: +852 2919-6441

**India:**  
India Medtronic Pvt. Ltd.  
Tel.: +(91)-80-22112245 / 32972359  
Mobile: +(91)-9611633007  
Patient Care Helpline: 1800 209 6777

**Indonesia:**  
Medtronic International Ltd.  
Tel.: +65 6436 5090 or +65 6436 5000

**Israel:**  
Medtronic Trading Ltd.  
Tel: +972-9-9724400  
Tel (product support - 8:00-17:00):  
+972-9-9724489  
Helpline: (weekends & holidays):  
1-800-611-888

**Italia:**  
Medtronic Italia S.p.A.  
Tel.: +39 02 24137 261  
Fax: +39 02 24138 210  
Servizio assistenza tecnica:  
N° verde: 800 60 11 22

**Japan:**  
Medtronic Japan Co. Ltd.  
Medtronic Japan  
24hr. Support Line: 0120-56-32-56  
日本: 日本メドトロニック株式会社  
24時間サポートライン: 0120-56-32-56

**Kazakhstan:**  
Medtronic BV in Kazakhstan  
Tel.: +7 727 321 13 30 (Aimaty)  
Круглосуточная линия поддержки:  
8 800 080 5001

**Kosovo:**  
Yess Pharma  
Tel.: +377 44 999 900  
Helpline: +37745888388

**Latin America:**  
Medtronic, Inc.  
Tel: 1(305) 500-9328  
Fax: 1(786) 709-4244

**Latvija:**  
RAL SIA  
Tel.: +371 67316372  
Helpline (9am to 6pm): +371 29611419

**Lithuania:**  
Monamedia UAB  
Tel.: +370 68405322  
Helpline: +370 68494254

**Macedonia:**  
Alkaloid Kons Dooel  
Tel.: +389 23204438

**Magyarország:**  
Medtronic Hungária Kft.  
Tel.: +36 1 889 0688

**Malaysia:**  
Medtronic International Ltd.  
Tel.: +603 7946 9000

**México:**  
Medtronic Servicios S. de R. L. de C.V.  
Tel. (México DF): +(11) 029 058  
Tel. (Interior): +01 800 000 7867  
Medtronic Directo 24/7 (from México DF): +(55) 36 869 787  
Medtronic Directo 24/7:  
+01 800 681 1845

**Middle East and North Africa:**  
Regional Office  
Tel.: +961-1-370 670

**Montenegro:**  
Glosarij d.o.o.  
Tel.: +382 20642495

**Nederland, Luxembourg:**  
Medtronic B.V.  
Tel.: +31 (0) 45-566-8291  
Gratis: 0800-3422338

**New Zealand:**  
Medica Pacifica  
Phone: 64 9 414 0318  
Free Phone: 0800 106 100

**Norge:**  
Medtronic Norge A/S  
Tel.: +47 67 10 32 00  
Fax: +47 67 10 32 10

**Philippines:**  
Medtronic International Ltd.  
Tel.: +65 6436 5090 or +65 6436 5000

**Россия:**  
ООО «Медтроник»  
Tel.: +7 495 580 73 77  
Круглосуточная линия поддержки:  
8 800 200 76 36

**Polska:**  
Medtronic Poland Sp. z o.o.  
Tel.: +48 22 465 6934

**Portugal:**  
Medtronic Portugal Lda  
Tel.: +351 21 7245100  
Fax: +351 21 7245199

**Puerto Rico:**  
Medtronic Puerto Rico  
Tel.: 787-753-5270

**Republic of Korea:**  
Medtronic Korea, Co., Ltd.  
Tel.: +82-2.3404.3600

**Romania:**

Medtronic Romania S.R.L  
Tel.: +40 372188017  
Helpline: +40 726677171

**Schweiz:**

Medtronic (Schweiz) AG  
Tel.: + 41 (0) 31 868 0160  
24-Stunden-Hotline: 0800 633333  
Fax Allgemein: +41 (0)318680199

**Serbia:**

Epsilon Research International d.o.o.  
Tel.: +381 113115554  
Medtronic Serbia D.o.o.  
Helpline: +381 112095900

**Singapore:**

Medtronic International Ltd.  
Tel.: +65 6436 5090 or +65 6436 5000

**Slovenija:**

Zaloker & Zaloker d.o.o.  
Tel.: +386 1 542 51 11  
24-urna tehnična pomoč:  
+386 51316560

**Slovenská republika:**

Medtronic Slovakia, s.r.o.  
Tel.: +421 26820 6942  
HelpLine: +421 26820 6986zilben

**Sri Lanka:**

Swiss Biogenics Ltd.  
Mobile: (+91)-9003077499  
or (+94)-777256760

**Suomi:**

Medtronic Finland Oy  
Tel.: +358 20 7281 200  
Help line: +358 400 100 313

**Sverige:**

Medtronic AB  
Tel.: +46 8 568 585 20  
Fax: +46 8 568 585 11

**Taiwan:**

Medtronic (Taiwan) Ltd.  
Tel.: 02-21836000  
Toll Free: +886-800-005285

**Thailand:**

Medtronic (Thailand) Ltd.  
Tel.: +662 232 7400

**Türkiye:**

Medtronic Medikal Teknoloji  
Ticaret Ltd. Sirketi.  
Tel.: +90 216 4694330

**Ukraine:**

ТОВ «Медтронік Україна»,  
Лінія цілодобової підтримки:  
Тел.: 0 800 508 300

**USA:**

Medtronic Diabetes Global  
Headquarters  
24-Hour Technical Support:  
+1-800-646-4633  
To order supplies: +1-800-843-6687

**United Kingdom:**

Medtronic Ltd.  
Tel.: +44 1923-205167

**Österreich:**

Medtronic Österreich GmbH  
Tel.: +43 (0) 1 240 44-0  
24 – Stunden – Hotline: 0820 820 190

