



Changing the MiniMed® Sure-T® Infusion Set with the MiniMed® 640G Insulin Pump



START HERE:

- Wash your hands. Press .
- Menu
Suspend Delivery
Audio Options
History
Reservoir & Tubing
Insulin Settings
Select **Reservoir & Tubing**.
- Reservoir & Tubing
New Reservoir
Fill Cannula
Select **New Reservoir**.
- New Reservoir
1. Remove infusion set from body.
2. Remove reservoir from pump.
Rewind
- Remove the infusion set you have been using by loosening the adhesive and pulling away from body.
- New Reservoir
1. Remove infusion set from body.
2. Remove reservoir from pump.
Rewind
Remove the used reservoir from the pump.

- New Reservoir
1. Remove infusion set from body.
2. Remove reservoir from pump.
Rewind
Select **Rewind**.
- New Reservoir
1. Fill reservoir.
2. Connect tubing to reservoir.
DO NOT CONNECT TO BODY.
Next

FILL RESERVOIR & CONNECT TO THE INFUSION SET TUBING
Follow the next steps to fill reservoir with insulin and connect to the infusion set tubing.

1. Remove from package. Make sure insulin vial is at room temperature to reduce the risk of air bubbles. (If using insulin cartridge proceed to step 6.)


If using insulin vial only. **If using insulin cartridge, skip to Step 6.**

2. Pull plunger down to the amount that you plan to fill with insulin.
3. Wipe vial with alcohol. Place vial on table. Firmly press the blue transfer guard onto vial.
4. Push and hold plunger down.
5. With your thumb still on the plunger, flip over so vial is on top. Release thumb and pull plunger down to fill with insulin.

Before you begin, please note that it is not recommended that you change your infusion set just prior to bedtime.

Changing your set during the day may reduce the risk of unexplained highs and no delivery alarms while sleeping. Please refer to the MiniMed® Sure-T® Infusion Set User Guide for more details.

To change your MiniMed® Sure-T® infusion set, you will need to organise the following supplies:



6. **If using a cartridge of insulin.** Hold the reservoir upright. Push plunger into the reservoir to expel any air. Turn plunger slightly anti-clockwise to loosen it. Firmly press insulin cartridge onto blue transfer guard. Hold the reservoir and cartridge at eye level and using a pencil, push down on the rubber stopper of the cartridge to fill the reservoir.
7. Tap the reservoir to move air bubbles to top of reservoir. Push plunger up to move air into vial.
8. If needed, pull plunger back down to amount of insulin needed for 2-3 days.
9. To avoid getting insulin on the top of the reservoir, turn vial over so it is upright. Hold transfer guard and turn reservoir counter-clockwise and remove from transfer guard.

IMPORTANT: If insulin or any liquid gets inside the tubing connector, it can temporarily block the vents that allow the pump to properly fill the infusion set. This may result in the delivery of too little or too much insulin, which could cause hypoglycaemia or hyperglycaemia.

CONNECT RESERVOIR TO INFUSION SET
You will place the reservoir connector onto the end of the infusion set to the filled reservoir.

1. Open infusion set packaging and connect the infusion set to the newly filled reservoir. Make sure both the top of the reservoir and the connector are dry before connecting them. Liquid can temporarily block the vents on the tubing connector.
2. Gently push connector onto reservoir. Turn clockwise until locked. You will hear a click.
3. If you see air bubbles, tap reservoir to move them to top. Push plunger just a bit to move them into tubing.
4. Twist plunger counter-clockwise to loosen and remove.



THE BACKLIGHT MAY HAVE TURNED OFF
Press any button to turn the screen back on.

Select **Load Reservoir** and unlock pump if necessary.

New Reservoir
1. Fill reservoir.
2. Connect tubing to reservoir.
DO NOT CONNECT TO BODY.
Next

PLACE RESERVOIR INTO PUMP
Now place the filled reservoir into the reservoir compartment of the pump.

New Reservoir
3. Place reservoir into pump and lock.
DO NOT CONNECT TO BODY.
Next

- Place reservoir into pump.
- Turn clockwise until you feel reservoir lock into place.
- New Reservoir**
3. Place reservoir into pump and lock.
DO NOT CONNECT TO BODY.
Next

Select **Next**.

LOAD RESERVOIR AND FILL TUBING
Follow these steps to load the reservoir and fill the tubing.

- Load Reservoir**
Select Load and hold until complete.
DO NOT CONNECT TO BODY.
Load Next
- Load Reservoir**
Complete
DO NOT CONNECT TO BODY.
Load Next
- Fill Tubing**
DO NOT CONNECT TO BODY.
Hold Fill until drops appear. Then select Next.
0.0
Fill Next
- Fill Tubing**
DO NOT CONNECT TO BODY.
Hold Fill until drops appear. Then select Next.
3.5
Fill Next

After you see drops, press and select **Next**.

Drops at end of tubing

INSERT INFUSION SET
Next, follow the steps to insert the infusion set into your body.

- Choose an insertion site from the shaded areas shown here. Wipe with alcohol or other antiseptic.
- Fill Cannula?**
1. Insert infusion set into body.
2. Select Fill to fill cannula or Done if not needed.
Fill Done
- Remove backing from the adhesive surrounding the needle.
- Remove needle guard.
- Insert infusion set using a 90 degree angle. Smooth down adhesive to secure set in place.
- Remove backing from adhesive at the tubing connection site. Place in a convenient location near insertion site. Ensure tubing between the two sites is not tight.

FILL CANNULA
When using a Sure-T Infusion Set you are not required to fill the canula. Please select "Done".

Select **Done**.

Fill Cannula?
1. Insert infusion set into body.
2. Select Fill to fill cannula or Done if not needed.
Fill Done

Your infusion set change is now complete!

