

# MINIMED™ 770G CGM TRAINING CHECKLIST

## PATIENT INFORMATION

NAME:	CERTIFIED PRODUCT TRAINER:	CLINIC NAME:	TRAINING DATE:
PUMP MODEL AND SERIAL NUMBER:	TRANSMITTER SERIAL NUMBER		

## KNOWLEDGE CHECK (understanding key terms and concepts)

<input type="checkbox"/> Completed MyLearning	<input type="checkbox"/> 24-hr Helpline 1800 777 808	<input type="checkbox"/> Difference between BG vs SG
<input type="checkbox"/> Auto CGM Getting Started Guide received	<input type="checkbox"/> Transmitter warranty (1 year)	<input type="checkbox"/> Focus on trends vs absolute SG values
<input type="checkbox"/> Medtronic apps installed	<input type="checkbox"/> How to charge transmitter (charging dock/green lights)	<input type="checkbox"/> Single, double and triple trend arrows
<input type="checkbox"/> CareLink™ Personal set up	<input type="checkbox"/> How to clean transmitter (green test plug)	<input type="checkbox"/> Fingertick confirmations required for treatment decisions
<input type="checkbox"/> Remove for X-ray, CT scan, MRI		

## PAIRING PUMP AND TRANSMITTER / SENSOR INSERTION

<input type="checkbox"/> Pair pump and transmitter	<input type="checkbox"/> Ensure thumb is on thumbprint when loading sensor
<input type="checkbox"/> Site selection, rotation and preparation	<input type="checkbox"/> Hold inserter against body after insertion (5-10secs)
<input type="checkbox"/> Sensor insertion technique	<input type="checkbox"/> Apply pressure to sensor adhesive

## TAPING (understanding key terms and concepts)

<input type="checkbox"/> Apply overtape to sensor (discuss options)	<input type="checkbox"/> Apply adhesive tab carefully onto transmitter
<input type="checkbox"/> Connect transmitter to sensor	<input type="checkbox"/> Apply overtape to transmitter

## STARTING THE SENSOR (understanding key terms and concepts)

<input type="checkbox"/> Start new sensor	<input type="checkbox"/> 3, 6, 12- and 24-hour graphs
<input type="checkbox"/> Sensor warm-up	<input type="checkbox"/> Sensor status screen
<input type="checkbox"/> Home screen and icons	

## PERSONALISE SETTINGS AND ALERTS (understanding key terms and concepts)

<b>SmartGuard™ low settings</b> <input type="checkbox"/> Time periods and limits <input type="checkbox"/> Alert before low <input type="checkbox"/> Alert on low	<b>Manual Mode only</b> <input type="checkbox"/> Suspend before low <input type="checkbox"/> Suspend on low <input type="checkbox"/> 2-hr max suspend <input type="checkbox"/> Auto vs manual resume	<b>SmartGuard™ high settings</b> <input type="checkbox"/> Time periods and limits <input type="checkbox"/> Alert before high <input type="checkbox"/> Time before high <input type="checkbox"/> Alert on high <input type="checkbox"/> Rise alert/Rise limit
<input type="checkbox"/> Snooze (duration of high and low)	<input type="checkbox"/> Alert Silence	<input type="checkbox"/> How to clear alerts

## CALIBRATION (understanding key terms and concepts)

<input type="checkbox"/> Importance of calibration and how to calibrate	<input type="checkbox"/> Avoid calibrating when 2-3 trend arrows are present
<input type="checkbox"/> Calibrate 3-4 times per day at stable times e.g. wake up, before main meals, before bed	<input type="checkbox"/> Calibrate immediately after BG test – do not wait
	<input type="checkbox"/> Calibration reminder (default is on)

## NOTES /FOLLOW UP APPOINTMENT SCHEDULE

Patient signature: \_\_\_\_\_ Date: \_\_\_\_\_

Trainer signature: \_\_\_\_\_ Date: \_\_\_\_\_