










What to do if SmartGuard™ Auto Mode is not ready









The **Auto Mode Readiness Table** shows what to do when the wait icon  or the question icon  appears by items on the Auto Mode Readiness screen. To open the Auto Mode Readiness screen, from the main menu, select **Status**, and then select **Auto Mode Readiness**.



Auto Mode Readiness Table

Line	Item	Instructions
1	Calibration required 	Perform a fingerstick and calibrate your sensor.
	BG required 	Perform a fingerstick and enter a new BG. Your BG must be within the 2.2–22.2 mmol/L range for your pump to enter Auto Mode.
	Wait to enter BG... 	Wait until the pump prompts you to enter a BG.
	Processing BG... 	Wait until the BG has processed.
2	Auto Mode turned off 	Turn on Auto Mode in the SmartGuard, Auto Mode screen.
3	Sensor not ready 	<p>a) Go to Utilities > Device Options to check if a transmitter ID has already been entered in your pump. For example, GT6133333F.</p> <p>If your pump does not have a transmitter ID entered, see <i>Pairing your pump and transmitter</i> in <i>GETTING STARTED WITH MINIMED™ 770G CONTINUOUS GLUCOSE MONITORING</i>.</p>

Quick Reference Guide for SmartGuard™ Auto Mode Readiness Screen

Line	Item	Instructions
		<p>b) Check your Home screen. If you see , move your pump and transmitter closer together. The pump will try to find the transmitter signal.</p> <p>If after 30 minutes the pump and transmitter are still not communicating, you will receive a Lost sensor signal alert. Check that the sensor is still inserted in the skin and the transmitter and sensor are still connected. Move your pump closer to your transmitter.</p>
		Turn on the sensor in the Utilities, Sensor Settings screen.
4		Wait until the bolus is complete or stop the bolus yourself before Auto Mode can activate.
5		If insulin delivery is suspended, Auto Mode cannot activate. Treat low BG if necessary as instructed by your healthcare professional.
6		When you turn on the Bolus Wizard feature for the first time, enter your Carb Ratio in the Edit Carb Ratio screen. You can also enter your Carb Ratio in the Bolus Estimate Setup screen, even if the Bolus Wizard feature is not turned on.
7		If a temp basal is currently active, you must wait until it has completed or cancel the temp basal yourself before Auto Mode can activate.
8		If active insulin is currently updating, it may take up to 5 hours to complete. You must wait until this amount is updated.
9		Auto Mode is gathering information on your insulin delivery history in order to personalize its automatic delivery of insulin.

Always read instructions for use