

Understanding the CareLink™ Adherence Report



Adherence (1 of 1)
08/04/2018 - 21/04/2018

one, patient

Glucose Measurements			Bolus Events				Fill Events					
BG Readings	Sensor Duration (h:mm)	Manual Boluses	Bolus Wizard Events	With Food	With Correction	Overridden	Rewind	Cannula Fills	Cannula Amount (U)	Tubing Fills	Tubing Amount (U)	
Sunday 08/04/2018	6	21:20	9	6	4		1	1	0.3	1	8.6	
Monday 09/04/2018	10	23:30	8	4	4							
Tuesday 10/04/2018	8	24:00	8	6	1							
Wednesday 11/04/2018	5	24:00	6	5	2							
Thursday 12/04/2018	3	24:00	8	6	1							
Friday 13/04/2018	8	24:00	9	5	4		1	1	0.3	1	10.2	
Saturday 14/04/2018	4	24:00	6	6	1							
Sunday 15/04/2018	6	21:10	8	6	3		1	1	0.3	1	3.2	
Monday 16/04/2018	10	24:00	11	7	4							
Tuesday 17/04/2018	8	24:00	12	9	4							
Wednesday 18/04/2018	4	24:00	9	7	3							
Thursday 19/04/2018	7	24:00	15	10	5		1	1	0.3	1	12	
Friday 20/04/2018	4	23:55	10	8	2							
Saturday 21/04/2018	1	8:05	3	3	1							
Summary	6.4/day	13d 02h 00m	0.0/day	9.2/day	72.1%	32.0%	0.0%	4	4	0.3U /fill	4	8.5U /fill


Partial day Suspend Suspend On Low Suspend Before Low Note: Partial days will not be included in summary averages. Days on which a time change occurred are considered to be partial days.

Medtronic

Adherence (1 of 1)
08/04/2018 - 21/04/2018




Take a look at your blood glucose readings every day. According to the American Diabetes Association, best practice is to check your blood glucose before each meal and before bedtime. That's about 4 BG readings per day. Do you show best practice? Do you meet the ADA's best practice recommendations?


1	Glucose Measurements		Bolus											Suspend Duration (h:mm)
	BG Readings	Sensor Duration (h:mm)	Ma Bol											
Sunday 08/04/2018	6	21:20												0:48
Monday 09/04/2018	10	23:30												0:15
Tuesday 10/04/2018	8	24:00												
Wednesday 11/04/2018	5	24:00												0:21
Thursday 12/04/2018	3	24:00												
Friday 13/04/2018	8	24:00												0:21
Saturday 14/04/2018	4	24:00		6	6	1								
Sunday 15/04/2018	6	21:10		8	6	3		1	1	0.3	1	3.2		0:05
Monday 16/04/2018	10	24:00		11	7	4								
Tuesday 17/04/2018	8	24:00		12	9	4								0:22
Wednesday 18/04/2018	4	24:00		9	7	3								0:38
Thursday 19/04/2018	7	24:00		15	10	5		1	1	0.3	1	12		0:08
Friday 20/04/2018	4	23:55		10	8	2								
 Saturday 21/04/2018	1	8:05		3	3	1								0:42
Summary	6.4/day	13d 02h 00m	0.0/day	9.2/day	72.1%	32.0%	0.0%	4	4	0.3U /fill	4	8.5U /fill		3h 34m

Take a look at your blood glucose readings every day. According to the American Diabetes Association, best practice is to check your blood glucose before each meal and before bedtime. That's about 4 BG readings per day. Do you show best practice? Do you meet the ADA's best practice recommendations?

Take a look at your Manual bolus (giving a set amount of insulin) vs Bolus Wizard™ Events. It is recommended you use the Bolus Wizard calculator when taking insulin for food and correcting high blood glucose readings. Using the bolus calculator benefits both you and your HCP in making easier adjustments in your pump settings. Manual boluses should be used only after discussions with your HCP.

Medtronic Adherence (1 of 1) one patient			08/04/2018			
			Glucose Measurements			
			BG Readings	Sense Duration (h:mm)		
Sunday	08/04/2018		6	21:20		
Monday	09/04/2018		10	23:30		
Tuesday	10/04/2018		8	24:00		
Wednesday	11/04/2018		5	24:00		
Thursday	12/04/2018		3	24:00		
Friday	13/04/2018		8	24:00		
Saturday	14/04/2018		4	24:00		
Sunday	15/04/2018		6	21:10		
Monday	16/04/2018		10	24:00		
Tuesday	17/04/2018		8	24:00		
Wednesday	18/04/2018		4	24:00		
Thursday	19/04/2018		7	24:00		
Friday	20/04/2018		4	23:50		
Saturday	21/04/2018		1	8:05		
Summary			6.4/day	13d 02h		

Be sure not to manually suspend your pump  for more than 1:00 hour unless advised by your HCP. Reasons to manually suspend your pump are commonly done for showering and water activities.

If you notice your pump suspended and you have a Suspend on Low , there was a time when your SG levels were low and your pump had to suspend your basal insulin. A SmartGuard™ Suspend event can last a maximum of 2 hours. If there are frequent occurrences of Threshold Suspend events, be sure to consult with your healthcare provider for possible changes to your basal rates.

		4	
		Insulin amount (U)	Suspend Duration (h:mm)
		8.6	0:48
			0:15
			0:21
		10.2	0:21
		3.2	0:05
			0:22
			0:38
		12	0:08
			0:42
		8.5U /fill	3h 34m

Understanding the CareLink™ Adherence Report

Medtronic

Adherence (1 of 1)
08/04/2018 - 21/04/2018

one, patient

	Glucose Measurements		Bolus Events					Fill Events					Suspend Duration (h:mm)
	BG Readings	Sensor Duration (h:mm)	Manual Boluses	Bolus Wizard Events	With Food	With Correction	Overridden	Rewind	Cannula Fills	Cannula Amount (U)	Tubing Fills	Tubing Amount (U)	
Sunday 08/04/2018	6	21:20		9	6	4		1	1	0.3	1	8.6	0:48
Monday 09/04/2018	10	22:00		9	6	4							0:15
Tuesday 10/04/2018	10	22:00		9	6	4							0:21
Wednesday 11/04/2018	10	22:00		9	6	4							0:21
Thursday 12/04/2018	10	22:00		9	6	4							0:05
Friday 13/04/2018	10	22:00		9	6	4							0:22
Saturday 14/04/2018	10	22:00		9	6	4							0:38
Sunday 15/04/2018	10	22:00		9	6	4							0:08
Monday 16/04/2018	10	22:00		9	6	4							0:42
Tuesday 17/04/2018	10	22:00		9	6	4							
Wednesday 18/04/2018	4	24:00		9	7	3							
Thursday 19/04/2018	7	24:00		15	10	5		1	1	0.3	1	12	
Friday 20/04/2018	4	23:55		10	8	2							
Saturday 21/04/2018	1	8:05		3	3	1							
Summary	6.4/day	13d 02h 00m	0.0/day	9.2/day	72.1%	32.0%	0.0%	4	4	0.3U /fill	4	8.5U /fill	3h 34m

5

Note how many times you Rewind your pump and how often. Remember you should be changing your reservoir and infusion set every 2-3 days for proper pump wear, so you should see a Rewind event occur every 2-3 days.

Medtronic

Medtronic Australasia Pty Ltd
2 Alma Road
Macquarie Park , NSW 2113
www.medtronic-diabetes.com.au

Safety Information: CareLink™ software

CareLink™ software is intended for use as a tool to help manage diabetes. The purpose of the software is to take information transmitted from insulin pumps, glucose meters and continuous glucose monitoring systems, and turn it into CareLink™ reports. The reports provide information that can be used to identify trends and track daily activities such as carbohydrates consumed, meal times, insulin delivery, and glucose readings. NOTE: CareLink™ report data is intended for use as an adjunct in the management of diabetes only and NOT intended to be relied upon by itself.

For a listing of indications, contraindications, precautions, warnings and potential adverse events please refer to the instructions for Use.

© 2018 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company.