













Diabetes product compatibility chart

Product	Transmitter	Sensor	Charging dock	Cleaning plug	Uploading device
MiniMed® 770G/780G 	Bluetooth Guardian Link 3†† 	Guardian Sensor 3# 	Grey Dock Charger takes 1 x AAA Battery. 	Green Test Plug 	Compatible Smartphone or Blue Adaptor≈ 
<p>i †† Transmitter has 'GL3' on front and Guardian Link 3 written underneath. The Bluetooth Guardian Link 3 will only work with Guardian Sensor 3.</p> <p># Sensors expire after 12 months and should be worn for a max of 7 days each.</p> <p>≈ For manual upload if automatic upload is not an option. As the 780G & 770G pumps automatically upload to CareLink, the new ACCU-CHEK® Guide Link Blood Glucose meter is not USB compatible. The MiniMed® 780G is indicated for use by patients age 7-80 years with Type 1 diabetes, whose total daily does of insulin is 8 units per day or more.</p>					













Product	Transmitter	Sensor	Charging dock	Cleaning plug	Uploading device
MiniMed® 670G 	Guardian Link 3** 	Guardian Sensor 3# 	Grey Dock Charger takes 1 x AAA Battery. 	Green Test Plug 	Contour Next Link 2.4 Only works with the 640G and 670G pumps. 
<p>i ** Transmitter has 'G' on front and Guardian Link 3 written underneath. The Guardian Link 3 will only work with Guardian Sensor 3.</p> <p># Sensors expire after 12 months and should be worn for a max of 7 days each.</p>					

Data is uploaded to: carelink.minimed.eu/app/login

Always follow the directions for use.

For detailed information regarding the instructions for use, indications, contraindications, warnings, precautions, and potential adverse events, please consult the device manual. For further information, contact your local Medtronic representative. ©2022 Medtronic. All Rights Reserved. Medtronic, MiniMed and Guardian are registered trademarks of Medtronic MiniMed, Inc. Please note: In contacting the Diabetes Toll Free, your personal and health information may be disclosed to an operator located outside Australia. 10677-022022

Diabetes product compatibility chart

Product	Transmitter	Sensor	Charging dock	Cleaning plug	Uploading device
MiniMed® 640G 	Guardian 2 Link* or Guardian Link 3** 	Enlite Sensor^ or Guardian Sensor 3# 	Grey Dock Charger takes 1 x AAA Battery. 	Grey test plug When using the Guardian Sensor 3, the cleaning plug will be Green not Grey. 	Contour Next Link 2.4 Only works with the 640G and 670G pumps. 
<div style="display: flex; align-items: center;"> <div style="font-size: 2em; margin-right: 10px;">i</div> <div> <p>* Transmitter has 'G' on front and Guardian 2 Link written underneath. The Guardian 2 Link will only work with Enlite Sensors.</p> <p>** Transmitter has 'G' on front and Guardian Link 3 written underneath. The Guardian Link 3 will only work with Guardian Sensor 3.</p> <p>^ Sensors expire after 6 months and should be worn for a max of 6 days each.</p> <p># Sensors expire after 12 months and should be worn for a max of 7 days each.</p> </div> </div>					
Guardian Connect Real Time CGM† 	Guardian Connect 3** 	Guardian Sensor 3# 	Grey Dock Charger takes 1 x AAA Battery. 	Green Test Plug 	iOS or Android Smartphone ⁰⁰ 
<div style="display: flex; align-items: center;"> <div style="font-size: 2em; margin-right: 10px;">i</div> <div> <p>† Phone not included.</p> <p>** Transmitter has 'GC' on front and Guardian Connect written underneath. Product codes 7821-Lxx (including letter "L") only work with Guardian Sensor 3.</p> <p># Sensors expire after 12 months and should be worn for a maximum of 7 days each.</p> <p>00 Upload automatically to cloud or manually upload via Carelink site.</p> </div> </div>					