

Empowering people on multiple daily injections to have **greater control of their diabetes.**^{1,2}

Guardian™ Connect

CONTINUOUS GLUCOSE MONITORING SYSTEM



1/2 THE WEIGHT
OF OTHER
AVAILABLE REAL-
TIME CGM^{3,4}

5.7g

**ACTUAL
SIZE**

35mm x 28mm



9mm

LOW PROFILE
SUPPORTS
AN ACTIVE
LIFESTYLE



WATERPROOF
UP TO 2.4m
FOR 30 MINUTE
PERIODS



**RE-
CHARGEABLE,**
REUSABLE,
LONG LIFE



**THE
AFFORDABLE
CGM⁵**

PER DAY^{*5}

\$9

*Based on 12-month supply.




FULLY REIMBURSED through the NDSS for eligible patients with T1D <21 years^{†6}



COST-EFFECTIVE for people transitioning from the NDSS^{5,7}

†Subsidised access to CGM sensors and transmitters is available through the NDSS to children and young people aged under 21 years, living with T1D, who meet specific criteria. Go to www.ndss.com.au/cgm for further information.



CONTROL FOR YOUR PATIENTS

CareLink®



SUPPORT FOR CARE PARTNERS AND HCPS

GUARDING THEIR PEACE OF MIND WITH PREDICTIVE ALERTS²

More advanced warning than
point-in-time alerts.

**THE ONLY CGM
AVAILABLE
WITH UP TO
1HR
PREDICTIVE
ALERTS**
of impending highs
or lows

**ENABLES
PEOPLE
TO AVOID
REBOUND**
events or
over-corrections



**ALLOWS
PARENTS OR
CARERS TO
RESPOND TO
CHILDREN.⁸**

Uses an algorithm based on the **rate of change** in glucose:

CURRENT
glucose
measure + **RATE**
of
change  = **PREDICTED**
time to
event

**For more information and downloads visit our dedicated HCP site:
hcp.medtronic-diabetes.com.au/guardian-connect**

References:

1. Liebl A *et al.* Continuous Glucose Monitoring: evidence and consensus statement for clinical use. *J Diabetes Sci Technol* 2013;7(2):500–19.
2. Pickup JC *et al.* Real-time continuous glucose monitoring in type 1 diabetes: a qualitative framework analysis of patient narratives. *Diabetes Care* 2015;38:544–50.
3. Guardian® Connect IFU.
4. Dexcom® G5 specifications.
5. www.medtronic-diabetes.com.au (accessed July 2017).
6. National Diabetes Services Scheme (NDSS). Continuous Glucose Monitoring Available at: <https://www.ndss.com.au/cgm> (accessed July 2017).
7. www.amsldiabetes.com.au/shop (accessed July 2017).
8. Seaquist ER *et al.* Hypoglycemia and diabetes: a report of a Workgroup of the American Diabetes Association and The Endocrine Society. *Diabetes Care* 2013;36(5):1384–95.

Medtronic

Medtronic Australasia Pty Ltd
5 Alma Road
Macquarie Park, NSW 2113
Australia
Tel: +61 2 9857 9000
Fax: +61 2 9889 5167
Toll Free: 1800 668 670

Medtronic New Zealand Ltd
Level 3 - Building 5, Central Park Corporate Centre
666 Great South Road
Penrose, Auckland 1051
New Zealand
Fax: +64 9 918 3742
Toll Free: 0800 377 807

E-Mail: australia.diabetes@medtronic.com
Facebook: www.facebook.com/MedtronicDiabetesAUS
Twitter: @DiabetesANZ
YouTube: Medtronic Diabetes ANZ

APPROVAL #4330-082017 © Medtronic 2017
All Rights Reserved. Printed in Australia 08/2017

medtronic.com.au
medtronic.co.nz