

FOCUS
ON YOUR
LIFE
**NOT
YOUR
LEVELS**



MiniMed™ 770G + dynamic algorithm.

The ONLY system that helps protect against
Lows AND Highs with Bluetooth connectivity.

Always read the instructions for use

Refers to SmartGuard™ Auto Mode. Some user interaction
required. Individual results may vary.

Medtronic

A PUMP THAT FOCUSES ON YOU, 24/7



*Components sold separately

EVERY 5 MINUTES

The MiniMed™ 770G insulin pump system* is engineered with advanced features that help keep you focused on your life.



PREDICTS

Monitors your glucose 24/7 and **predicts** whether you're going high or low.



NO ASSUMPTIONS

No pre-set basal rate means **no assumptions** about insulin requirements.³



SELF-ADJUSTS

Smooths out highs and lows, with SmartGuard™ tech that increases or decreases insulin delivery, **every 5 minutes.**




ADAPTS

Every day at midnight, the dynamic algorithm recalculates and resets the total insulin dose. It **adapts** to how your body works.²

MiniMed™ 770G THE ONLY SYSTEM THAT AUTOMATICALLY ADJUSTS LEVELS FOR YOU

Worry less with **SmartGuard™ tech** that self-adjusts, taking on more tasks of managing diabetes for you - so you can keep your day on track.

SMARTGUARD™ AUTO MODE



HELPS PREVENT HIGHS

If you're going high, it gives you more insulin

HELPS PREVENT LOWS

If you're going low, it gives you less insulin




"The ONLY system that helps protect against Lows AND Highs, increasing Time in Range"^{1,2}



Refers to SmartGuard™ Auto Mode. Some user interaction required. Individual results may vary.

STAY IN THE KNOW ON THE GO

Don't get stalled during busy or active days—the MiniMed™ 770G system keeps you informed with real-time glucose trends on your phone, so you can glance and go.

-  SmartGuard™ ON
-  Low predicted, insulin suspended
-  Glucose in range



MiniMed™ Mobile App



**KNOW IF YOU GO
HIGH OR LOW**

Get notified on your phone if you're going high or low. At any time, glance at the app to instantly know your real-time glucose levels. Auto-Synch to Carelink reports to support in-office or remote doctor visits

CareLink™ Connect App



**GIVE LOVED ONES
PEACE OF MIND**

Securely share pump and CGM info with up to five care partners.



Current CGM reading

Insulin delivery

CGM trend over time



Check for compatible devices at www.medtronic-diabetes.com.au/minimed-770g-supported-devices



WIRELESS UPGRADES

Wireless Bluetooth Software enables future updates and upgrades to the latest technology from Medtronic**

**Subject to eligibility

MiniMed™ 770G Insulin Pump

Works as the hub for the system. The pump communicates with the other system components to calculate insulin adjustments.

Infusion Set

Allows for continuous and discreet delivery of insulin.

Continuous Glucose Monitoring (CGM)

Guardian™ Sensor 3 measures glucose levels every 5 minutes, sending info to the pump.

Smartphone App

Get notified on your phone if you are going high or low. Easily track your pump and CGM info with the MiniMed™ Mobile app.

Blood Glucose Meter





Accurate test results from the Accu-Chek® Guide Link blood glucose meter are sent wirelessly to the pump for quick sensor calibrations.

Estimated availability in the first quarter of 2021.



LEAVE ASSUMPTIONS BEHIND.

Move forward with the dynamic algorithm.

	 VEO™	 MiniMed™ 640G	 MiniMed™ 670G Auto Mode	 MiniMed™ 770G Auto Mode
	2010	2015	2019	2020
▶ Self-adjusting basal rate	✗ Pre-set	✗ Pre-set	✓ Calculated every 5 minutes	✓ Calculated every 5 minutes
▶ Self-adjusting total daily dose/ Max basal rate	✗ Pre-set by HCP	✗ Pre-set by HCP	✓ Automatically re-calculated every midnight	✓ Automatically re-calculated every midnight
▶ Self-adjusting insulin sensitivity factor	✗ Pre-set by HCP	✗ Pre-set by HCP	✓ Automatically re-calculated every midnight	✓ Automatically re-calculated every midnight
▶ Helps protect against highs and lows	✗ Lows only	✗ Lows only	✓ Highs and lows	✓ Highs and lows
▶ Wireless Software Upgrades	✗	✗	✗	✓ Via Bluetooth
▶ Glucose Readings on Smart Phones	✗	✗	✗	✓ Via Bluetooth
▶ Automatic Syncing with Carelink	✗	✗	✗	✓ Via compatible Smart Phone

v HCP=healthcare professional.

Based on your patient's needs

ALWAYS READ THE INSTRUCTIONS FOR USE (IFU)

For detailed information regarding indications, contraindications, warnings, precautions, and potential adverse effects, please consult the IFU. Always consult a Health Care Professional before making treatment decisions, and to see if this product is suitable for you.

Medtronic

Medtronic Australasia Pty Ltd
2 Alma Road Macquarie Park,
NSW 2113 Australia

Tel: +61 2 9857 9000
Fax: +61 2 9889 5167
Toll Free: 1800 668 670

medtronic-diabetes.com.au

Email: australia.diabetes@medtronic.com
Facebook: www.facebook.com/MedtronicDiabetesAUS
YouTube: [Medtronic Diabetes ANZ](https://www.youtube.com/MedtronicDiabetesANZ)

APPROVAL: #8444-092020 © Medtronic 2020
All rights Reserved. Printed in Australia 06/2020

References:

1. Bergenstal RM. et al. JAMA. 2016;316(13):1407–1408.
2. Garg SK et al. Diabetes Technol Ther. 2017;19(3):155–163.
3. Protocol for hybrid closed loop technology. MiniMed™ 670G system. Medical education.