




Home Screen

Press

- Suspend Delivery
- Audio Options
- History
- Reservoir & Tubing
- Insulin Settings
- Sensor Settings
- Event Markers
- Reminders
- Utilities



Medtronic Diabetes 24-hour toll free helpline (Australian landlines): **1800 777 808**

Medtronic Australasia Pty Ltd, 97 Waterloo Road, North Ryde NSW 2113 Australia. *Only available when sensor feature is on.
 Please note: In contacting the Diabetes Toll Free, your personal and health information may be disclosed to an operator located outside Australia.
 ©2014 Medtronic Australasia Pty Ltd. All rights reserved. EX14-102014



Home Screen

Press

Suspend Delivery
Audio Options
History
Reservoir & Tubing
Insulin Settings
Sensor Settings
Event Markers
Reminders
Utilities

Previously Reservoir + Set

Previously Sensor

All reminders found in one place

Reservoir & Tubing
New Reservoir
Fill Cannula

Settings
Sensor
Sensor Connections*
Alert Silence*
High Settings*
Low Settings*
Calibrate Sensor*
Auto Calibration*

Reminders
Personal
Bolus BG Check
Missed Meal Bolus
Low Reservoir
Set Change
Calibration*

Previously Suspend, now first item on Menu

Suspend Delivery

All History found in one place

History
Summary
Daily History
Alarm History

Previously Capture Events

Event Markers
BG
Injections
Food
Exercise
Other

Much like current Utilities menu

Utilities
Airplane Mode
Display Options
Time/Date
Remote Bolus
Block
Self Test
Carb Unit
Manage Settings
Sensor Demo
Device Options
Keypad Lock
Language

Previously Alert + Type

Audio Options
Audio
Audio & Vibrate
Vibrate
Volume

Settings that effect insulin delivery are found here

Insulin Settings
Bolus Wizard Setup
Basal Pattern Setup
Preset Temp Setup
Preset Bolus Setup
Dual/Square Setup
Bolus Increment
Max Basal/Bolus
Easy Bolus
Auto Suspend
Bolus Speed



Medtronic Diabetes 24-hour toll free helpline (Australian landlines): **1800 777 808**